


PASSIVE MATTRESSES

LINET

Portfolio



For comfort and safety

Quality mattresses from LINET for all your needs




Wide selection according to the patient's condition




A broad spectrum of foam types and mattress designs for the various needs of hospital wards and nursing homes. Choose the right surface for your patient/client on the following pages.

Safety and comfort




The mattresses are designed to improve pressure redistribution through the mechanisms of immersion and envelopment. They also offer high levels of lying and sleeping comfort.

High-quality materials and production

















Only the most high-grade materials are used in a production process that meets international standards. Therefore, the mattresses remain fully functional for a long lifetime. Quality management system certified according to DIN EN ISO 9001 and 13485. Environmental management system certified according to DIN EN ISO 14001.

Infection control




LINET mattresses are protected by a waterproof and vapour permeable cover. The cover is antibacterial equipped with long durability and stability. The surface and sides are made of proven PERLASTIC® Silber** according to the guidelines for bed hygiene in clinics and care.

Three segments of passive mattresses

	Recommended areas of use	Features	Foam	Cover
PREMIUM	<ul style="list-style-type: none">— Acute care— ICU and critical areas (ViskoMatt 50) 	<p>A combination of the highest quality foams, especially for high-risk areas.</p> <ul style="list-style-type: none">— Efficient pressure injury prophylaxis— Excellent comfort for lying and sleeping— Special cold-foam in the gluteal area for extra long durability – even with high patient weight loads 	<p>Viscoelastic foam is often referred to as “memory foam”. It reacts to body heat and body weight. This allows the body to sink into the mattress slightly, providing optimal support and adjustment to the body’s contours. The resulting pressure relief is the best choice for pressure injury prevention.</p> 	<p>PERLASTIC® Silber protective cover</p> <ul style="list-style-type: none">— High infection prophylaxis— Handles for an optimized handling— Hidden extra strong zipper— Bi-elastic and dimensionally stable, long durability and stability— Extra strong PU coating for permanent protection against bacteria and viruses— Breathable— Antimicrobial equipped— Latex free  <p>A perfect fit of the cover: no wrinkles for excellent prophylaxis of pressure injuries.</p>
ADVANCED	<ul style="list-style-type: none">— Acute care— Nursing care for the elderly— Long-term care 	<p>Combination of cold and standard PU foam for better pressure relief designed for patients at medium risk of pressure injuries.</p> <ul style="list-style-type: none">— With edge zone for self-mobilization of the patient.— Good pressure relief and high comfort— Safety for patients 	<p>A cold foam is superior quality material due to its open pores. It ensures high air permeability and pleasant sleeping conditions. It quickly returns to its original shape without permanent deformation.</p> 	  <p>waterproof antibacterial</p>
VALUE	<ul style="list-style-type: none">— Acute care— Nursing care for the elderly— Long-term care 	<p>Entry level monoblock passive mattresses designed for low risk patients.</p> <ul style="list-style-type: none">— Profiling for improved immersion and envelopment— High comfort for sleeping— Durability 	<p>A standard PU foam cut to alter the deformation properties of the foam. These types of modifications break the surface tension of the foam increasing its ability to conform to the body shape.</p> 	<p>The mattress cover is a close fit to the foam core inside, protecting both the mattress and the patient’s skin. The cover includes a number of technical features that supports this dual function.</p> <ul style="list-style-type: none">— Water resistant— Multi-stretched— Vapour permeable— Fungal and microbacterial resistance   <p>waterproof antibacterial</p>


Passive Mattresses Portfolio

			VALUE			ADVANCED	PREMIUM				
			<div>EffectaCare 10</div> <div>— Monoblock foam mattress with wave profiling on both sides which makes it flippable inside the cover.</div> <div><div>sewn seams110 kgdouble-sided</div></div>	<div>EffectaCare 20</div> <div>— Monoblock foam mattress. The core is cut transversely on top in the sacral area. Bottom profiling allows easier mattress positioning.</div> <div><div>sewn seams150 kg</div></div>	<div>PrimaCare 10</div> <div>— One layer mattress with transverse and longitudinal profiling at the top part. Fully profiled for additional contouring.</div> <div><div>sewn seams200 kg</div></div>	<div>MediMatt 30</div> <div>— Two lying surfaces with high elastic foam for high comfort. Lying surface with a polyether-cold foam and support layer from polyether foam.</div> <div><div>sewn seams150 kg</div></div>	<div>ViskoMatt 10</div> <div>— Entry-level model for pressure injury prevention with good pressure relief. Viscoelastic foam surface with a polyether support layer.</div> <div><div>welded seams150 kg</div></div>	<div>ViskoMatt 30</div> <div>— Very good pressure relief. Highly elastic and durable with an antimicrobial viscoelastic polyurethane foam surface and a polyether cold foam support layer.</div> <div><div>welded seams165 kg</div></div>	<div>ViskoMatt 50</div> <div>— Excellent pressure relief. Highly elastic and durable with an antimicrobial viscoelastic polyurethane foam surface and a polyether cold foam support layer.</div> <div><div>welded seams190 kg</div></div>	<div>ViskoMatt 10B</div> <div>— Mattress with higher maximum patient weight. Excellent pressure relief. Highly elastic and durable, combined of 3 layers: an antimicrobial viscoelastic polyurethane foam layer, and 2 layers with polyether foam.</div> <div><div>sewn seams235 kg</div></div>	<div>ViskoMatt 30B</div> <div>— Mattress for bariatric patients. Excellent pressure relief. Highly elastic and durable, combined of 3 layers: an antimicrobial viscoelastic polyurethane foam layer, polyether cold foam layer and polyether layer.</div> <div><div>sewn seams285 kg</div></div>
			Monoblock foam			2 layer foam	Viscoelastic foam				
Foam Layer	Upper					5 cm cold foam 50 kg/m³	5 cm viscoelastic foam 52 kg/m³	5 cm viscoelastic foam 85 kg/m³	7 cm viscoelastic foam 85 kg/m³	5 cm viscoelastic foam 85 kg/m³	6 cm viscoelastic foam 85 kg/m³
	Middle	14 cm monoblock Polyether foam 28 kg/m³	14 cm monoblock Polyether foam 33 kg/m³	14 cm monoblock Polyether foam 41 kg/m³	×	×	×	×	×	5 cm Polyether foam 41 kg/m³	6 cm cold foam 45 kg/m³
	Lower				9 cm Polyether foam 41 kg/m³	9 cm cold foam 50 kg/m³ & 33 kg/m³	9 cm cold foam 80 kg/m³ & 33 kg/m³	6 cm Polyether foam 46 kg/m³	8 cm Polyether foam 46 kg/m³		
Reinforced Edge		×	×	×	✓	✓	✓	✓	✓	✓	✓
Zipper		180°	180°	180°	180°	180°	180°	180°	180°	180°	180°
Cut Corners		×	×	✓	✓	✓	✓	✓	✓	✓	✓
Ergoframe		×	×	×	×	×	×	×	×	×	×
Profiled Top Surface		Convolutad	Partial	Full	×	×	×	×	×	×	×
Max. Patient Weight		110 kg	150 kg	200 kg	150 kg	150 kg	150 kg	165 kg	190 kg	235 kg	285 kg
Min. Patient Weight		40 kg	40 kg	40 kg	30 kg	30 kg	30 kg	30 kg	30 kg	30 kg	30 kg
2 Transport Handles		×	×	×	✓	✓	✓	✓	✓	4 transport handles	4 transport handles
Evacuation Cover (option)		×	×	×	✓	✓	✓	✓	✓	✓	×
Foam/Cover Fire Protection		Crib 5	Crib 5	Crib 5	Crib 5	Crib 5	Crib 5	Crib 5	Crib 5	Crib 5	Crib 5



Foam connection without adhesive seams


- Uniform distribution of the contact pressure without hard seams
- Highest pressure relief



Edge Zone Reinforcement

Strong 10 cm reinforcement at the peripheral zone:

- Easy self-mobilization of the patient
- Excellent fall protection



Labelling

UDI compliant labelling


- According to the requirements of the european MDR
- Cover and core labelled for easy traceability

ViskoMatt Premium Mattresses

ViskoMatt 50

The upper layer of the ViskoMatt 50 is made of viscoelastic foam to ensure a high level of comfort and pressure injury prevention. The mattress also features


a reinforced pelvic area to prevent pressure injuries when the backrest is raised.



Visko Matt 30

ViskoMatt 30 offers optimum pressure relief of 35% on average. The high-quality foams ensure high stability and


excellent durability. In addition to this, local pressure points are avoided.



ViskoMatt 10

ViskoMatt 10 has viscoelastic foam to provide optimal pressure distribution for all parts of the body that are in

contact with the mattress, considerably reducing the risk of pressure injuries.



Why Viscoelastic = Memory Foam

- The material “remembers” the burden and after getting up the contours of the body are still visible.
- The core of viscoelastic mattresses adapt to each individual body form depending on temperature and weight.
- The pressure is spread equally and the body is ideally relieved = mode of action for pressure injury prophylaxis.



Why cold foam

- Cold foam (HR, High resilience) is created by technologically modifying standard polyurethane foam. It is produced at lower temperatures to ensure that it has superior properties.
- It has an open-cell structure that makes it very breathable
 - It is durable, flexible, and dimensionally stable with high point elasticity



Quality defined by density



- The higher the density of foam, the:
- greater the restoring force and pressure release
 - longer the durability of the mattress


Density specifies how many kilograms of foam there are per cubic meter.
Example: density 60 = 60 kg/m³

PERLASTIC® Silber protective cover

For MediMatt 30, ViskoMatt 10, 10B, 30, 30B, 50. Total weight 230 g/m².

The trickot sock between the cover and the core for a long mattress lifetime

The cover consists of two parts: outside PU-coating and inside polyester-textile




Extra stable spiral-zipper

Individual handles



PU-Coating

- Extremely durable and abrasion-resistant. Triple coating with 90g/m².
- Bacteria and virus proof (Liquid-tight according to ISO 1420 > 200 cm)
- Tight against fluids and constant against urine, blood and ointments
- Antimicrobial, contains zinc ions which prevent the development of mould
- Breathable (Water vapor transmission at ASTM 96-66 550 g/m² in 24 hour by 38 °C)




With the tricot sock, which is between cover and core, the cover can easily be changed. It reduces friction and shearing forces and increases the mattress lifetime.

Textile ensures

the stability and the PU coating the function. The textile is an especially fine woven product with narrow mesh that provides high durability and reduces the elasticity to decrease wrinkles. High grammage 140 g/m².

Option: Evacuation Cover Pro

Best risk management: simple and fast evacuation – combined with the highest hygiene standards.



Extra-wide, all-round lip protects the carrying system from contamination.



Long carrying handles at the head and foot end can be unfolded quickly and enable an uncomplicated one-man evacuation.

Infection Control

- Full protective cover: lying surface and sides made of proven PERLASTIC® Silber
- Tucked away straps: Protection against germs in hard-to-reach places

Safe Evacuation

- 2 optimally positioned straps for securing the patient
- Wide belt buckles for optimal adjustment
- Suitable for patients weighing up to 250 kg (valid for the evacuation straps)

Mattresses Dimensions (cm)

EffectaCare 10	200 × 86 × 14	190 × 86 × 14	ViskoMatt 30	200 × 86,5 × 14	208 × 96 × 14
	200 × 90 × 14	190 × 90 × 14		200 × 90 × 14	208 × 96 × 14
EffectaCare 20	200 × 86 × 14	190 × 90 × 14		190 × 86,5 × 14	208 × 88 × 14
	200 × 90 × 14	142 × 70 × 10		190 × 90 × 14	200 × 114 × 14
	190 × 86 × 14	137 × 70 × 10		200 × 100 × 14	200 × 79 × 14
	200 × 100 × 16			220 × 86,5 × 14	190 × 79 × 14
PrimaCare 10	200 × 86 × 14	190 × 86 × 14	ViskoMatt 50	200 × 86,5 × 16	208 × 96 × 16
	200 × 90 × 14	190 × 90 × 14		200 × 90 × 16	208 × 96 × 16
MediMatt 30	200 × 86,5 × 14	208 × 96 × 14		190 × 86,5 × 16	208 × 88 × 16
	200 × 90 × 14	208 × 96 × 14		190 × 90 × 16	220 × 86,5 × 16
	190 × 86,5 × 14	208 × 88 × 14		200 × 100 × 16	200 × 114 × 16
	190 × 90 × 14	200 × 114 × 14	ViskoMatt 10B	200 × 86,5 × 16	
	200 × 100 × 14	200 × 79 × 14		208 × 86,5 × 16	
ViskoMatt 10	220 × 86,5 × 14	190 × 79 × 14		200 × 90 × 16	
	200 × 86,5 × 14	208 × 96 × 14		200 × 100 × 16	
	200 × 90 × 14	208 × 88 × 14		200 × 120 × 16	
	190 × 86,5 × 14	142 × 70 × 10		200 × 114 × 16	
	190 × 90 × 14	137 × 70 × 10		205 × 114 × 16	
	208 × 96 × 14	200 × 114 × 14	ViskoMatt 30B	200 × 86,5 × 20	
	200 × 100 × 14	200 × 79 × 14		208 × 86,5 × 20	
	220 × 86,5 × 14	190 × 79 × 14		200 × 90 × 20	
				200 × 100 × 20	



Members of LINET Group

Bigla Care AG | Bernstrasse 3 | CH-3421 Lyssach
T +41 31 700 92 00 | info@bigla-care.ch | www.bigla-care.ch